GradiPlasma LA Low Lupus Anticoagulant Quality Control Plasma, Low Level

for use with LA Screen/LA Confirm (DRVVT)

For *in vitro* diagnostic use.

INTENDED USE
GradiPlasma LA Low is a low positive control plasma for use in Lupus Anticoagulant clotting test assays, specifically LA Screen (DRVVT) and LA Confirm (DRVVT) from Life Diagnostics.

SUMMARY AND EXPLANATION

"Lupus Anticoagulant" (LA) was a term used in the early 1970s to describe a condition whereby non-specific circulatory coagulation inhibitors were found to be present particularly in cases of systemic lupus erythematosis (SLE)¹. Currently LA is recognized as being associated with a variety of thrombotic syndromes².

REAGENTS

For *in vitro* diagnostic use.

Caution: Treat as potentially infectious.

This product has been collected from lupus anticoagulant patients under FDA approved conditions. All human plasma prepared for lyophilization have been found negative when tested for antibodies to HIV-1/2 and HCV. Non-reactive for hepatitis B surface antigen (HBsAg) by FDA required tests.

Pack contents
10 vials lyophilized platelet poor human plasma from LA positive and normal donors, buffers and stabilizers, sodium azide.

When disposing of sodium azide, always flush with large volumes of water to avoid the possibility of an explosive residue forming in metal plumbing.

Additional materials required but not provided
Purified water, USP or equivalent
Micropipet, 1.0 ml

Storage conditions
Store at 2-8 °C
The lyophilized plasma should remain stable at least until the expiry date stated on the vial when stored at 2-8 °C Following reconstitution, plasma should be used within 4 hours at +20 °C open vial) or 8 hours at 2-8 °C (closed vial). Reconstituted plasma may be frozen for 1 week at -20 °C.

TEST PROCEDURE

NB: All standard precautionary practices (hygiene, spill procedures, disposal of biological materials etc) must be observed.

- Reconstitute with the volume of purified water stated on the vial label.
- Mix gently, do not shake.
- Allow to stand at room temperature for at least 15 minutes before use.
- Prepare Lupus Anticoagulant testing reagents according to manufacturers specifications.
- Test undiluted QC plasmas at the same time as patient samples.
- Record results of clotting times and/or ratios on Quality Control charts as required.

PERFORMANCE CHARACTERISTICS
GradiPlasma LA Low should give a ratio between 1.4 and 1.7 using a Dade Behring BCS instrument with Life Diagnostics' LA Screen/Confirm test system. Other reagent/instrument combinations may yield different results.

A new range based on mean + 2SD should be established with each lot of controls and reagents in each test system.


Within run and total precision were calculated for the GradiPlasma LA Low by testing twice daily over 20 days using an IL MLA800 instrument. Within run and total precision for both volumes of control were less than 3% with LA Screen and LA Confirm reagents.

Validation studies for the GradiPlasma LA Low have been performed using IL ACL300, IL MLA 800 and Dade Behring BCT instruments.

AVAILABILITY OF GRADIPLASMA LA PRODUCTS
LAHP-1 10 x 1ml vials GradiPlasma LA High
LALP-1 10 x 1ml vials GradiPlasma LA Low

BIBLIOGRAPHY

- Feinstein D.I., Rapaport, S.P. Acquired inhibitors of blood coagulation. Progress in haemostasis and thrombosis, Vol 1. Ed Spael T., New York, Grune and Stratton, 75-95, 1972
- Schleider M.A. et al. A clinical study of the lupus anticoagulant. Blood 48: 499-509, 1976

	
--	---

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