## LA Screen / LA Confirm

### Simplified Dilute Russell’s Viper Venom Test (DRWT) for detection of Lupus Anticoagulants

#### INTENDED USE

LA Screen and LA Confirm are simplified DRVVT reagents for detection of Lupus Anticoagulants (LA) in one-stage clotting tests.

#### LA Screen

Simplified DRVVT reagent to screen for the presence of Lupus Anticoagulants.

#### LA Confirm

Phospholipid-rich DRVV reagent for the specific correction of Lupus Anticoagulants.

#### SUMMARY AND EXPLANATION

LAs are autoantibodies against negatively charged phospholipids or complexes of phospholipids with either beta-2-glycoprotein 1 or clotting factors such as prothrombin. They occur in various clinical conditions, especially autoimmune diseases<sup>1</sup> and are now considered to be a significant risk factor in patients with otherwise unexplained thrombosis and are often present in women who have recurrent fetal loss<sup>2</sup>. LA have traditionally been detected using phospholipid responsive clotting tests, such as the activated partial thromboplastin time (APTT), kaolin clotting time (KCT) and DRVVT<sup>2</sup> where they have an anticoagulant effect.

The DRVVT time first became popular following the publication by Thiagarajan et al. in 1986<sup>3</sup>. Life Diagnostics has simplified and standardized this method.<sup>4</sup> According to Petri<sup>5</sup> et al, thrombosis in SLE patients is more closely linked with LA detectable by DRVVT than with anticardiolipin antibodies detected by enzyme linked immunosorbent assays (ELISA)<sup>5</sup>. This observation was recently extended by Galli and Bevers<sup>6</sup> who showed that the LA subtype with most effect on DRVVT tests is beta-2-glycoprotein 1 dependent and different to the prothrombin-requiring LA subtype having more prolonging effect on KCT tests.

#### Test principle

- Russell's viper venom present in LA Screen initiate plasma clotting by directly activating factor X. LA antibodies prolong the LA Screen clotting time.
- LA Confirm reagent is similar to LA Screen but contains a high phospholipid concentration. The extra phospholipid counteracts the LA antibody and largely corrects the clot time<sup>3</sup>.
- DRVV tests "bypass" factor VII of the extrinsic pathway and the contact and antithaemophilic factors of the intrinsic pathway. Therefore LA Screen is more specific for LA than APTTs as they are not affected by contact factor abnormalities or by factor VIII deficiencies or antibodies<sup>5</sup>. There have been a number of tests developed based on phospholipid correction<sup>4</sup> however none have been as convenient to use as LA Confirm subsequent to LA Screen.
- Mixing tests may be useful to exclude factor II, V and X deficiencies, which may prolong LA Screen and LA Confirm results. Mixing normal plasma with test plasma replenishes any factors deficient in the test plasma. If the mixing test is still prolonged, it indicates that an inhibitor (such as LA) is present in the test plasma.

#### REAGENT

##### Composition

LA Screen and LA Confirm contain Russell's viper venom, phospholipids, antiheparin agents, calcium, buffers, stabilizers, sodium azide and dyes

#### Warnings and precautions

For in-vitro diagnostic use only.

When disposing of azide, always flush with large volumes of water to avoid the possibility of an explosive residue forming in metal plumbing.

#### Preparation for use

Reconstitute with the volume of distilled water stated on the vial. Mix well by inversion to ensure complete re-suspension of the lyophilized material.

#### Storage instructions

The lyophilized reagents are stable at least until the expiry date stated on the vial label when stored at 2-8°C.

Following reconstitution reagents can be stored for 8 hrs at 37°C, 24 hours at 20-25°C, 48 hrs at 2-8°C and 1 month at -20°C.

#### Indications of instability / deterioration

If there is no evidence of vacuum when the vial is opened and/or the reagent does not appear dry it should be returned to Life Diagnostics or your local distributor.

#### SPECIMEN COLLECTION AND PREPARATION

Collect and process blood in accordance with NCCLS Standard H21-A3: Collection, Transport and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline-3rd Edition (1998). Separated plasma should be stored at 2-8°C and tested within 4 hours of collection. If the samples are to be frozen before testing, the plasma must be double spun to remove platelets (to below 10 x 10<sup>9</sup>/L)<sup>10</sup> prior to freezing as these can shorten the LA Screen clotting time.

#### TEST PROCEDURE

- Pre-warm a slight excess of LA Screen or LA Confirm reagent at 37±1°C in a reagent reservoir, allowing for 200µl per test.
- Dispense 200µl of test plasma into a glass test tube and warm for 1 minute at 37°C
- Add 200µl of pre-warmed reagent to the plasma and time from the moment of reagent addition to a clotting end point.
- Repeat for duplicate test values and report the average of these as the result.

#### Automated method

Protocols for most automated clotting machines are available on request. LA Screen tests should be carried out using equal volumes of test plasma and reagent as in most thrombin time protocols.

However, observation or acquisition times should be extended to 120 seconds.

#### Materials provided

Each pack of LA Screen/ LA Confirm  
 Order code LASD-10 contains: 10 x lyoph. for 2 ml  
 Order code LASD-25 contains: 10 x lyoph. for 5 ml  
 Order code LACD-5 contains: 10 x lyoph. for 1 ml  
 Order code LACD-10 contains: 10 x lyoph. for 2 ml

#### Materials required but not provided

10mm x 75mm glass test tubes  
 Waterbath at constant temperature of 37 +/- 1°C  
 Stopwatch  
 200µl precision pipette  
 Purified water, USP or equivalent  
 Platelet depleted normal plasma  
 1, 2 or 5ml pipette.  
 Normal and Abnormal plasma for Quality Control

#### Quality control

Each laboratory should determine its own acceptable control values and normal range. Gradiplasma LA High (code LAHP) and Gradiplasma LA Low (code LALP) are abnormal control plasmas manufactured by Life Diagnostics for use in quality control of LA Screen and LA Confirm assays.  
 Quality control plasma must be tested at the same time as patient samples.

#### RESULTS


If the LA Screen clotting time is within the normal range no further testing for LA may be necessary. If the LA Screen clotting time result is more than 20% longer than pooled platelet depleted normal plasma the result should be considered abnormal and investigated further.

(See Figure 1)

The final result is best expressed as a ratio of the clotting times of LA Screen divided by LA Confirm.

LA Ratio =  $\frac{\text{LA Screen clotting time}}{\text{LA Confirm clotting time}}$

FIGURE 1

#### Mixing tests

If results are borderline, mixing studies may be used to correct for hidden defects in the sample and clarify the presence of LA. These tests should be carried out on a 50:50 mixture of test plasma and platelet depleted normal plasma using the standard test procedure.

TABLE 1: Combination of mixing and confirmatory tests

LA Screen Clotting time		LA Confirm Clotting Time		Interpretation
Patient Plasma	Mix- Patient + Normal	Patient Plasma	Mix- Patient + Normal	
N	N	N	N	LA not detected
ABN	ABN	N	N	LA probably present
ABN	N	ABN	N	Possible factor deficiency/OAT exclude by further investigation
ABN	ABN	ABN	N	Possible factor deficiency exclude by further investigation
ABN	ABN	ABN	ABN	Exclude other inhibitor by further investigation

#### LIMITATIONS OF THE PROCEDURE

Samples containing clots and those with abnormal hematocrit should be discarded. Jaundiced, lipemic and hemolyzed specimens should be tested by manual techniques as some photoelectric instruments give false results.

Commercially available normal quality control plasmas with unspecified levels of citrate and platelets are not recommended for use in mixing studies.  
 For comparative studies LA Screen and LA Confirm tests should be performed at the same time.

LA assays based on different properties appear to be more or less sensitive to certain subgroups of LAs. Therefore at least two screening assays, based on different properties, should be performed before the possibility of LA is excluded.<sup>10</sup>

Heparin levels up to 1 unit/ml have no effect due to the presence of a neutralizing agent in both LA Screen and LA Confirm.

#### EXPECTED VALUES

Normal ranges for LA Screen of 31-44 seconds and for LA Confirm of 30-38 were obtained using a manual tilt-tube method with samples collected from 26 healthy individuals aged 18 to 55 years.

The ratio of LA Screen / LA Confirm was in the range 0.8-1.2. These results should be used as a guide only.

Each laboratory should determine its own normal range for both manual and automated methods to overcome differences in sample collection and instrumentation.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

Intra-laboratory and inter-laboratory precision studies were performed using both manual tilt-tube and automated methods. These studies gave coefficient of variation less than 3.5% on normal plasmas and less than 5% for abnormal samples.

Specificity studies were performed on known plasma samples.


A positive ratio for LA Screen / LA Confirm of greater than 1.2 was found in known\*

LA plasma -	90% (26/29)
Heparinised plasma-	12% (1/8)
OAT patient plasma-	0% (0/7)
Factor deficient plasmas-	0% (0/5)
Normal plasma-	2% (1/60)

\*N.B. "known" LA identification with APTT + KCT mixing studies.

#### BIBLIOGRAPHY

- Feinstein DJ. Lupus anticoagulant, thrombosis and fetal loss. N.Eng.J.Med. 1985, 313,1248-1350
- Triplett DA, Brandt JT. Laboratory identification of the lupus anticoagulant. Br.J. Haematol. 1991,73,139-142
- Exner T, Triplett DA, Taberner D, Machin SJ. Guidelines for testing and revised criteria for lupus anticoagulants. Thromb. Haemostas. 1991,65,320-322
- Rauch J, Tannenbaum M, Janoff AS. Distinguishing plasma lupus anticoagulants from antifactor antibodies using hexagonal phase II phospholipids. Thromb. Haemostas. 1989,62,892-896
- Love PE, Santoro SA. Antiphospholipid antibodies: Anticardiolipin and the lupus anticoagulant in systemic lupus erythematosus (SLE) and in non-SLE disorders. Ann. Intern. Med. 1990,112,692-698
- Thiagarajan P, Pengo V, Shapiro SS. The use of the dilute Russell's viper venom time for the diagnosis of lupus anticoagulants. Blood. 1986,68,869-874
- Exner T, Papadopoulos G, Kouits J. Use of a simplified dilute Russell's viper venom time (DRVVT) confirms heterogeneity among "lupus anticoagulants". Bl. Coag. Fibrinol. 1990,1,259-266
- Petri M, Rheinschmidt. Metal. The frequency of lupus anticoagulant in systemic lupus erythematosus. Ann. Int. Med. 1987,106,524-531
- Galli M, Finazzi G, Bevers EM, Barbui T. Kaolin clotting time and dilute Russell's viper venom time distinguish between prothrombin dependent and beta-2-glycoprotein 1 dependent antiphospholipid antibodies. Blood. 1995,86,617-623
- Brandt JT, Triplett DA, Alving B, Scharrer I. Criteria for the Diagnosis of Lupus Anticoagulants: An Update. Thromb.Haemostas. 1995, 74 (4), 1185-90

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#### BIBLIOGRAPHIE

- Feinstein DJ. Lupus anticoagulant, thrombosis and fetal loss. N.Eng.J.Med. 1985, 313,1248-1350
- Triplett DA, Brandt JT. Laboratory identification of the lupus anticoagulant. Br.J. Haematol. 1991,73,139-142
- Exner T, Triplett DA, Taberner D, Machin SJ. Guidelines for testing and revised criteria for lupus anticoagulants. Thromb. Haemostas. 1991,65,320-322
- Rauch J, Tannenbaum M, Janoff AS. Distinguishing plasma lupus anticoagulants from antifactor antibodies using hexagonal phase II phospholipids. Thromb. Haemostas. 1989,62,892-896
- Love PE, Santoro SA. Antiphospholipid antibodies: Anticardiolipin and the lupus anticoagulant in systemic lupus erythematosus (SLE) and in non-SLE disorders. Ann. Intern. Med. 1990,112,682-698
- Thiagarajan P, Pengo V, Shapiro SS. The use of the dilute Russell's viper venom time for the diagnosis of lupus anticoagulants. Blood, 1986,68,869-874
- Exner T, Papadopoulos G, Kouits J. Use of a simplified dilute Russell's viper venom time (DRVVT) confirms heter